



Anaesthetic breathing systems

Dr Ami Kotecha




Introduction

- 3 main functions
 - Delivery of anaesthetic gases and vapours
 - Oxygenation of the patient
 - Removal of carbon dioxide
- Several different classifications exist



Properties of the ideal breathing system

- Simple and safe to use
- Delivers intended inspired gas mixture
- Permits spontaneous, manual & controlled in all ages
- Efficient, requiring low flow rates
- Protects patients from barotrauma
- Sturdy, compact and lightweight
- Easy removal of waste exhaled gases
- Easy to maintain with minimal running costs




Old classification

Open	No boundary and no dead space	Oxygen tubing near patient
Semi-open	Partial boundary between airway and atmosphere	Schimmelbusch mask
Semi-closed	Fully bounded. Prevents entry of atmospheric air but vents excess fresh gas	Mapleson systems
Closed	No venting of excess gas	Circle systems at low flows

Schimmelbusch mask







Old classification

Open	No boundary and no dead space	Oxygen tubing near patient
Semi-open	Partial boundary between airway and atmosphere	Schimmelbusch mask
Semi-closed	Fully bounded. Prevents entry of atmospheric air but vents excess fresh gas	Mapleson systems
Closed	No venting of excess gas	Circle systems at low flows

Mapleson systems







Old classification

Open	No boundary and no dead space	Oxygen tubing near patient
Semi-open	Partial boundary between airway and atmosphere	Schimmelbusch mask
Semi-closed	Fully bounded. Prevents entry of atmospheric air but vents excess fresh gas	Mapleson systems
Closed	No venting of excess gas	Circle systems at low flows


Circle system





Modern classification



Systems with potential for rebreathing	Bidirectional flow within the system	Mapleson classification
Non-rebreathing systems	Unidirectional flow with valves Carbon dioxide absorption	Drawover systems and resuscitation bags Circle systems and 'to and fro' systems




Components of a breathing system

- Adjustable pressure-limiting valve
- Reservoir bag
- Tubing

Adjustable pressure-limiting valve




- Spill valve, pop-off valve, expiratory valve, relief valve
- Designed to vent gas during positive pressure
- Pressure of less than 0.1 kPa activates the valve when open




Components

- 3 ports: inlet, patient and exhaust port-latter can be open to atmosphere or connected to the scavenging system
- Lightweight disc sits on a knife-edge seating- held in place by a spring
- Tension in the spring and therefore the valve's opening pressure is controlled by the valve dial



Mechanism of action

- One-way, adjustable, spring-loaded valve
- Valve allows gases to escape when the pressure in the breathing system exceeds the valve's opening pressure
- During spontaneous ventilation, the patient generates a positive pressure during expiration, causing the valve to open
- During positive pressure ventilation, a controlled leak is produced in inspiration by adjusting the valve dial, allowing control of the patient's airway pressure




Problems in practice and safety features

- Malfunction of scavenging system can cause excessive negative pressure which can lead to the APL valve being open throughout respiration
- The patient can be exposed to excessive positive pressure if the valve is shut during assisted ventilation- safety mechanism actuated at 6 kPa in some designs
- Water vapour in exhaled gas may condense on the valve- disc made of hydrophobic material which prevents water condensing

Reservoir bag




- Antistatic rubber or plastic
- Ellipsoid in shape
- Standard adult size is 2 l (range from 0.5 to 6 l)



Mechanism of action

- Accommodates fresh gas flow during expiration acting as a reservoir available for the following inspiration
- Acts as a monitor of patient's ventilatory pattern during spontaneous breathing and also a very inaccurate guide to tidal volume
- Used to assist or control ventilation




Problems in practice and safety features

- Because of its compliance, the reservoir bag can accommodate rises in pressure in the breathing system
- Limits the pressure in the system to about 4 kPa
- Size of the bag depends on breathing system and patient
- If too small may not provide a large enough reservoir for big tidal volumes
- If too large, does not act as respiratory monitor

Tubing




- Corrugated or smooth
- Different lengths depending on system being used
- Allow humidification of inspired air
- Parallel and coaxial arrangements available




Mapleson classification

- Mapleson classified breathing systems in 1954 into 5 types (later modified to include the Rees modification)
- Efficiency: fresh gas flow required to prevent rebreathing
- Rebreathing:
 - Rebreathing expired air including carbon dioxide
 - Recirculation of expired gas with carbon dioxide removed and oxygen added


Mapleson classification



Mapleson A




- Corrugated rubber or plastic tubing: 110- 130 cm in length
- Reservoir bag at machine end
- APL valve at the patient end




Mechanism of action

- As patient exhales, dead space gas is channelled towards reservoir bag which is continuously filled with fresh gas
- Pressure builds, opening APL valve and expelling alveolar gases
- The patient inspires, getting a mixture of fresh gas and dead space gases
- Efficient for spontaneous breathing: FGF required to prevent rebreathing = patient's minute volume (70 ml/kg/min)
- Not efficient for controlled ventilation: FGF required to prevent rebreathing = x2-3 minute volume

Mechanism of action







Problems in practice and safety features


- Not suitable for children less than 30 kg due to increased dead space at patient end
- Heaviness of APL valve at patient end puts a lot of drag on the connections especially if connected to scavenging system

Mapleson A- Lack system




- Coaxial modification of Magill Mapleson A
- 1.8 m length
- FGF through outside tube, exhaled gases through inner tube
- Inner tube wide in diameter (14 mm) to reduce resistance to expiration
- Reservoir bag at machine end
- APL valve at machine end eliminating drag at the patient end

Mapleson B and C





- Components
 - Reservoir bag
 - APL valve at patient end
 - FGF proximal to APL valve
- Not efficient for spontaneous ventilation (FGF of x1.5-2 MV required to prevent rebreathing)
- During controlled ventilation, B system more efficient as tubing acts as reservoir (FGF > 50% still required to prevent rebreathing)

Mapleson D- Bain circuit



- Coaxial and parallel versions available (180 cm)
- FGF through inner tube, exhaled gas through outer tube
- Reservoir bag at machine end
- APL valve at machine end





Mechanism of action

- During spontaneous ventilation, patient's exhaled gases are channelled back to the bag and mixed with FGF
- Pressure build up opens the APL valve allows venting of this mixture of gases
- FGF required to prevent rebreathing is about 150-200 ml/kg/min
- Inefficient for spontaneous ventilation
- More efficient for controlled ventilation
- FGF of 70 ml/kg/min

Mechanism of action





Figure 4



Mode of action of Mapleson D breathing system during spontaneous ventilation.




Problems in practice and safety features

- Internal tube can kink preventing FGF to the patient
- Internal tube can become disconnected at the machine causing a huge increase in dead space resulting in hypoxaemia and hypercapnia
- Movement of the reservoir bag is not always an indication that fresh gas is being delivered to the patient

Mapleson E and F




- Valveless breathing system used for children upto 30 kg
- Suitable for spontaneous and controlled ventilation
- Components
 - T shaped tubing with 3 ports
 - FGF delivered to one port
 - 2nd port goes to patient
 - 3rd port leads to reservoir tubing (Jackson-Rees modification)




Mechanism of action

- Requires a FGF of x2-3 minute volume to prevent rebreathing (minimum flow 4 l/min)
- Double-ended bag acts as a monitor for ventilation
- Bag can be used for assisted or controlled ventilation and CPAP during spontaneous ventilation
- Volume of the reservoir tubing determines degree of rebreathing (too large a tube) or entrainment of room air (too small a tube)
- Volume of reservoir bag should be about patient's tidal volume



Problems in practice and safety features

- No APL valve so scavenging can be a problem
- High flows required to prevent rebreathing
- Double-ended bag allows CPAP (particularly useful in children under 6 years of age as they have a low FRC)



Topics to look up

- Checking circuits- Bain's
- Circle circuits
 - Components of a circle circuit
 - Draw a circle circuit
- Humphrey ADE breathing system
- Drawover systems